

DAFTAR LAMPIRAN

Lampiran 1. Perhitungan Jumlah Dehidrator	50
Lampiran 2. Alat dan Bahan yang digunakan pada penelitian	51
Lampiran 3. Dokumentasi Penelitian	52
Lampiran 4. Kromatogram MEFA Minyak Jarak	54
Lampiran 5. Kromatogram MEFA DCO P ₂ O ₅ 1 % w/w; T = 150°C; t = 2 jam <i>Variasi Jumlah dehidrator</i>	55
Lampiran 6. Kromatogram MEFA DCO P ₂ O ₅ 1,5 % w/w; T = 150°C; t = 2 jam <i>Variasi Jumlah dehidrator</i>	56
Lampiran 7. Kromatogram MEFA DCO P ₂ O ₅ 2 % w/w; T = 150°C; t = 2 jam <i>Variasi Jumlah dehidrator</i>	57
Lampiran 8. Kromatogram MEFA DCO P ₂ O ₅ 2,5 % w/w; T = 150°C; t = 2 jam <i>Variasi Jumlah dehidrator</i>	58
Lampiran 9. Kromatogram MEFA DCO P ₂ O ₅ 3 % w/w; T = 150°C; t = 2 jam <i>Variasi Jumlah dehidrator</i>	59
Lampiran 10. Kromatogram MEFA DCO P ₂ O ₅ 3 % w/w; T = 100°C; t = 2 jam <i>Variasi Suhu</i>	60
Lampiran 11. Kromatogram MEFA DCO P ₂ O ₅ 3 % w/w; T = 125°C; t = 2 jam <i>Variasi Suhu</i>	61
Lampiran 12. Kromatogram MEFA DCO P ₂ O ₅ 3 % w/w; T = 150°C; t = 2 jam <i>Variasi Suhu</i>	62
Lampiran 13. Kromatogram MEFA DCO P ₂ O ₅ 3 % w/w; T = 175°C; t = 2 jam <i>Variasi Suhu</i>	63
Lampiran 14. Kromatogram MEFA DCO P ₂ O ₅ 3 % w/w; T = 200°C; t = 2 jam <i>Variasi Suhu</i>	64
Lampiran 15. Kromatogram MEFA DCO P ₂ O ₅ 3 % w/w; T = 200°C; t = 0,5 jam <i>Variasi Waktu</i>	65
Lampiran 16. Kromatogram MEFA DCO P ₂ O ₅ 3 % w/w; T = 200°C; t = 1 jam <i>Variasi Waktu</i>	66

Lampiran 17. Kromatogram MEFA DCO P ₂ O ₅ 3 % w/w; T = 200°C; t = 1,5 jam <i>Variasi Waktu</i>	67
Lampiran 18. Kromatogram MEFA DCO P ₂ O ₅ 3 % w/w; T = 200°C; t = 2 jam <i>Variasi Waktu</i>	68
Lampiran 19. Kromatogram MEFA DCO P ₂ O ₅ 3 % w/w; T = 200°C; t = 2,5 jam <i>Variasi Waktu</i>	69
Lampiran 20. Kromatogram MEFA DCO P ₂ O ₅ 3 % w/w; T = 200°C; t = 3 jam <i>Variasi Waktu</i>	70
Lampiran 21. Kromatogram MEFA DCO P ₂ O ₅ 3 % w/w; T = 200°C; t = 1 jam <i>Variasi Waktu</i>	71
Lampiran 22. Kromatogram MEFA DCO P ₂ O ₅ 3 % w/w; T = 200°C; t = 2 jam <i>Variasi Waktu</i>	72
Lampiran 23. Kromatogram MEFA DCO P ₂ O ₅ 3 % w/w; T = 200°C; t = 3 jam <i>Variasi Waktu</i>	73
Lampiran 24. Kromatogram MEFA DCO P ₂ O ₅ 3 % w/w; T = 200°C; t = 4 jam <i>Variasi Waktu</i>	74